



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C3orf22-chromosome 3 open reading frame 22 Gene

Reference: PTXBC032025

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

uncharacterized protein C3orf22; chromosome 3 open reading frame 22

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 22
GENE ID	152065
ORIGIN SPECIES	Human
ACCESSIONS	BC032025
VECTOR	pDONR223
NCBI SYMBOL	C3orf22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagactcagtgctgctgaagaagctcaccagagtaagaagtgtaggagatccagccggcaggagaattttgccaaagaattctcgatagagtgctgtgctgacagagcccgaccctgagccctgagccctggagggtgcacaaacgactcgaacacgggtgcagctgcccctgcaga
agaggtggttgccaaagagtgctaccgagctccgagggctcggggctgcagatttatacatcagctggtgctctgcggccacatcagacacatcaccttccatgctgcaacctttgggaactcaagttctgtagtgcgccctctccagacaactcgcttctgctgtgcacc
cgacacatgctgctgctgcgccccagacagcagcgccgagcgtctgcaggggctccctctcgt