



Reference: PTXBC032100

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCA57; TID1; hTID-1; dnaJ homolog subfamily A member 3, mitochondrial; DnaJ (Hsp40) homolog, subfamily A, member 3; dnaJ protein Tid-1; hepatocellular carcinoma-associated antigen 57; tumorous imaginal discs protein Tid56 homolog; DnaJ heat shock protein		

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily A, member 3
GENE ID	9093
ORIGIN SPECIES	Human
ACCESSIONS	BC032100
VECTOR	pDONR223
NCBI SYMBOL	DNAJA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1