

JAGN1-jagunal homolog 1 (Drosophila) Gene

Reference: PTXBC032101

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GL009; SCN6; protein jagunal homolog 1; jagunal homolog 1		

For research use only.

GENE DESCRIPTION	jagunal homolog 1 (Drosophila)
GENE ID	84522
ORIGIN SPECIES	Human
ACCESSIONS	BC032101
VECTOR	pDONR223
NCBI SYMBOL	JAGN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtctcgagcagggccgagcggccgacccgagcgagcgttccagcaccgggagcggcgtcgccatgcactaccagatgagtgtagccctcaagtatgaaatcaagaagctgatctacgtacatctggctcatatggctgctggtgctaagatgagcgtgggacacct
gaggctcttgcacatgatcaggtggccatgccctatcagtggaataccggtatttgcagcatttgcctctctcttggccttctctcttccccgcaacaacattagctacctggctcctcatgatcagcagtgaggactctttccatcgtccactcattatggcagcatggagatgttc
cctgctgcacagcagctctaccgcatggcaaggcctaccgtttcctcttgggttttctgcggttccatcatgtacctggtgtgtgtgtggcagtgcaagtgcatgcctggcagttgtactacagcaagaagctcctagactcttgggtcaccagcacacaggagaagaagcataaatga