



Reference: PTXBC032108

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VPS52, GARP complex subunit; ARE1; SAC2; SACM2L; dJ1033B10.5; vacuolar protein sorting-associated protein 52 homolog; SAC2 suppressor of actin mutations 2-like protein; vacuolar protein sorting 52 homolog		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 52 homolog (S. cerevisiae)
GENE ID	6293
ORIGIN SPECIES	Human
ACCESSIONS	BC032108
VECTOR	pDONR223
NCBI SYMBOL	VPS52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr