

C20orf11-chromosome 20 open reading frame 11 Gene

Reference: PTXBC032120

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein C20orf11; C20orf11; TWA1; glucose-induced degradation protein 8 homolog; two hybrid associated protein No. 1 with RanBPM; two hybrid-associated protein 1 with RanBPM; GID complex subunit 8 homolog		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 11
GENE ID	54994
ORIGIN SPECIES	Human
ACCESSIONS	BC032120
VECTOR	pDONR223
NCBI SYMBOL	C20orf11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttatgcagaaaaaccgatgaaatcacgaaagatgagtggaagctcaataactgcatgtccagagacgagacatgaaccgctcatgaactacgtgtcacagagggtttaagggaagcagcggagaagttcgaatggaatctggaatcgaactagtttgatct
ggaaacacttgatgaacgaatcaagatccgggagatgatactgaaaggtcagattcaggaggccatcgcttgatcaacagcctccaccagagctctggacacaaaccggatcttacttccatttgcagcaacagcattgatcgagctgatccgagcgggagacagagcgccg
gctggagtttcacagactcagctggcggagcgggagagagccgagagtgccctcacagagatggagcgtaccctggcactgctggcctttgacagtcccaggagctgccccttcggagacctctccacaccatgcagaggcagaaggtgtgagtgaaagtaaccaagctgt
gctagattatgaaatcgagtgcaacccaaactggcaaaattactgaaactacttgggctcagaacgagctggaccagaagaaagtaaaatccccaaatgacagacctcagcaagggtgtgattgaggagcccaagtag