

DDX39-DEAD (Asp-Glu-Ala-Asp) box polypeptide 39 Gene

Reference: PTXBC032128

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DDX39; BAT1; BAT1L; DDXL; URH49; ATP-dependent RNA helicase DDX39A; DEAD (Asp-Glu-Ala-Asp) box polypeptide 39; DEAD (Asp-Glu-Ala-Asp) box polypeptide 39A; DEAD box protein 39; DEAD-box helicase 39A; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 39; UAP56-r		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 39
GENE ID	10212
ORIGIN SPECIES	Human
ACCESSIONS	BC032128
VECTOR	pDONR223
NCBI SYMBOL	DDX39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaacaggatgtgaaaaacgatcttttgattacgatgaagaggaagagccccaggctcctcaagagagcacaccagctccccctaagaagacatcaaggatcctacgtttccatccacagctctggctccgggactttctgctgaagccggagctcctgcgggcatcgtg
gactgtggctttgagcatcctctgaggtccagcatgagtgattccccaggccatcctgggcatggagctcctgtgcccaggccaagtcgggatgggcaagacagcgggtcttcgtgctggccacctacagcagattgagcctgtcaacggacaggtgacggctcgtgcatgtgccac
acgaggagctggtcctccagatcagcaaggaatagagcgtttccaagtacatgccagcgtcaaggtgtctgttcttcgggtgtctctccatcaagaagatgaagaagtgttgaagaagaactgtcccatgtcgtgggtgggacccccggccgcatcctggcgctcgtgcgga
ataggagcttcagcctaagaatgtgaagcactttgtgctggacgagtgacaaagatgctggagcagctggacatgcggcgggatgtgcaggagatcttcgcctgacaccacacgagaagcagtcagatgttcagcgccacctgagcaaggacatccggcctgtgtgcaggaa
gttcagtcaggatcccatggaggtgtttgtggagcagagagaccaagctcacgctgcacggcctgcagcagtagtactacgtcaaaactcaagacagtgagaagaaccgcaagctcttgaatctgtgctggtgagtttaaccagcctgtcacgctatcagcagttcaaggattccagcgg
cggatcctggtggccacaaatctgtttggccggggatggacatcagcagtgatcaacatcgtctttaa