

TIMM10-translocase of inner mitochondrial membrane 10 homolog (yeast) Gene

Reference: PTXBC032133

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TIM10; TIM10A; TIMM10A; mitochondrial import inner membrane translocase subunit Tim10; translocase of inner mitochondrial membrane 10 homolog; translocase of inner mitochondrial membrane 10		

For research use only.

GENE DESCRIPTION	translocase of inner mitochondrial membrane 10 homolog (yeast)
GENE ID	26519
ORIGIN SPECIES	Human
ACCESSIONS	BC032133
VECTOR	pDONR223
NCBI SYMBOL	TIMM10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcctctcaggcccaacagctggctgcggagctggaggtggagatgatggccgatgtacaacagaatgaccagtgcctgccacgggaagtgtgtgcctcctcactacaaggaagcagagctctccaagggcgagctgtgtgcctggaccgatgtgtcttaagtacctggac
atccatgagcggatgggcaaaaagttgacagagttgtctatgcaggatgaagagctgatgaagaggggtgcagcagagctctggcctcatga