

## C9orf23-chromosome 9 open reading frame 23 Gene

Reference: PTXBC032136

|                |                                                                                                                                                                                                                                   |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | alba-like protein C9orf23; C9orf23; bA296L22.5; ribonuclease P protein subunit p25-like protein; RNase P protein subunit-like p25; ribonuclease P/MRP 25kDa subunit-like; rpp25-like protein; ribonuclease P/MRP subunit p25 like |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 9 open reading frame 23                                                                |
| <b>GENE ID</b>          | 138716                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC032136                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C9orf23                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggagcactaccgaaagctggctctgtagagctccagcgccttcccaatgccagctacctctgatacccttgagatgcgggtccgagatggcagcaaaattcgcaacctgctggggtggctctgggtcggttgaggcgccagctgctggcatgtagttctcaggtctg  
gcagggtgcaggaaagctgtcagctgcctgagattgtcaagcggcggtccagcctgcaccagctcaccaagctactttccttcagactgaggacagctgggtccagcctcactgacacaggcctagacccctcacagtgcgcccgcctatgtcctgcagtggtgctg  
gtcagccgggacccctggaccccaatgagtggtgtaccaaccccgaggacccccctggcctgggttccatgccagctccagctgtggcctctgtccgaagaaggctcgagaccccgatcgtga