

FGFR1OP2-FGFR1 oncogene partner 2 Gene

Reference: PTXBC032143

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC123-like; WIT3.0; FGFR1 oncogene partner 2; fibroblast growth factor receptor 1 oncogene partner 2; wound inducible transcript 3.0		

For research use only.

GENE DESCRIPTION	FGFR1 oncogene partner 2
GENE ID	26127
ORIGIN SPECIES	Human
ACCESSIONS	BC032143
VECTOR	pDONR223
NCBI SYMBOL	FGFR1OP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttgcacaattgagaagcacttgcgcgctaaagctctgttgaagattaagagatcatgacgatgcagcagaatctctgattgagcaaaccacagctctcaacaagcgagtagaagccatgaacagatcaggagaattcaagaacttaagaagtcgcgagacatcgcc
cacgggtccaggttagttatgggaatccagcaagaaaacagacaaatcagagagttgcaacaagaaaacaagaattacgtacatctctggaagaacatcagtcggccttggaacttataatgagcaagtaccgagaacaaatgtttagattgctaattggctagcaaaaagatgatccgggt
ataataatgaagttaaaagagcagcactccaagattgacatggtacatcgtacaagtcggaaggattctctctgatgcacatcgcacacatcctgaagcacctcaacatggactggagagaaggcacttgggaagcaaatcagaatgtacactaa