

TAC3-tachykinin 3 Gene

Reference: PTXBC032145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HH10; NKB; NKNB; PRO1155; ZNEUROK1; tachykinin-3; gamma tachykinin 3; neurokinin b; neuromedin K; preprotachykinin B; tachykinin 3		

For research use only.

GENE DESCRIPTION	tachykinin 3
GENE ID	6866
ORIGIN SPECIES	Human
ACCESSIONS	BC032145
VECTOR	pDONR223
NCBI SYMBOL	TAC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggatcatgctgtatcacagccatcctggccttcagcctagctcagagcttgggctgtctgtaaggagccacaggagggtgttcttggcggggccgcagcaagaggatccagatctctaccagctgtccagagactcttcaaaagccactcatctctgagggttgc
tcaaagccctgagccaggctagcacagatcctaaggatcaacatctcccagaaaacgtgacatgcatgacttcttgggacttatggcgaagaggagcgtccagccagactctctacggatgtgaatcaagagaacgtccccagcttggcatcctcaagtatccccgagagcaga
atag