

TNNI2-troponin I type 2 (skeletal, fast) Gene

Reference: PTXBC032148

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AMCD2B; DA2B; FSSV; fsTnI; troponin I, fast skeletal muscle; troponin I fast twitch 2; troponin I type 2 (skeletal, fast); troponin I, fast-twitch isoform; troponin I, fast-twitch skeletal muscle isoform; troponin I, skeletal, fast; troponin I2, fast skel		

For research use only.

GENE DESCRIPTION	troponin I type 2 (skeletal, fast)
GENE ID	7136
ORIGIN SPECIES	Human
ACCESSIONS	BC032148
VECTOR	pDONR223
NCBI SYMBOL	TNNI2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagatgaggagaagcggacaggccatcagggccgcagggcagcactgaagagcgtgatgctgcagatagcggccacggagctggagaaggaggagagccgctgaggcagagaagcagaactacctggcggagcactgcccgcctgcataccgggtcca
tgcctgaagtcaggagctctgcaaacagctgcacgccaagatcgcggctgaagaggagaagtacgacatggagggtgagggtgcagaagaccagcaaggagctggaggacatgaaccagaagctatttgatctcggggcaagttcaagcggcccccactgcgagggtgctg
catgtcggccgatgccatgtcaaggccctgctgggctgaagcacaagggtgtgcatggacctgagggccaacctgaagcaggtcaagaaggaggacacagagaaggagcggacctgcgagacgtgggtgactggagaagaacatcgaggagaagtctggcatggaggcc
ggaagaagatgtttgagtcgagtcctag