

MGC40069-hypothetical protein MGC40069 Gene

Reference: PTXBC032242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC40069
GENE ID	348035
ORIGIN SPECIES	Human
ACCESSIONS	BC032242
VECTOR	pDONR223
NCBI SYMBOL	MGC40069
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctggagcttatccactgtctgggatacatgttctctgagaactgccagagccagtcagtgacccagcctgacatccacatcactgtctctgaaggagcctcactggagttgagatgaactatctctatggggcaacaccttatctcttgatggaaagaactgttgaagaagct
ttttactctgtgttgcctgaagccttgagggtagcaagtagtctggagaaaaaggagaaggagatgaatcttccagctgctacttggtcacgctataatgtattagaggcacactgtctctacctttgataagatggctaacctcaggtgactcactgttatcagctcagccacattgccca
caggaactgtga