

CCDC88C-coiled-coil domain containing 88C Gene

Reference: PTXBC032316

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAPLE; HKRP2; KIAA1509; SCA40; protein Daple; Dvl-associating protein with a high frequency of leucine residues; hook-related protein 2; spinocerebellar ataxia 40; coiled-coil domain containing 88C		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 88C
GENE ID	440193
ORIGIN SPECIES	Human
ACCESSIONS	BC032316
VECTOR	pDONR223
NCBI SYMBOL	CCDC88C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggggcgccgaatgggaaagcgaacagactgggtgaaaactttggcccggttggaagcggcagcaggacacacgactatgtacatggatttagtggacggcatcttttgaacaaattatgctgcaaatagatcccaggcccaaatcaacgcataataagcagtcaca
atgatgtgaaccttcgcattcagaatttgacctcttggtgagaacattaagacactactaccagatagaccttttccggtaa