

KIAA1704-KIAA1704 Gene

Reference: PTXBC032333

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1704; AD029; LSR7; bA245H20.2; GPALPP motifs-containing protein 1; lipopolysaccharide specific response-7 protein; lipopolysaccharide-specific response protein 7; GPALPP motifs containing 1		

For research use only.

GENE DESCRIPTION	KIAA1704
GENE ID	55425
ORIGIN SPECIES	Human
ACCESSIONS	BC032333
VECTOR	pDONR223
NCBI SYMBOL	KIAA1704
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaagagacgtgatggaccggccctgccgcccgttcaaggcccgaggacgaggaagcgggaccgagccctgttcaggaccagctctgccccctaatataaaagcagtagttcagattcatcagacagcgaagacagtagttctgtacgaagaag
gaaatcaagaatctgaagaagatgacagtggtcccaactgcaagaaaacagaggaaaatcaggatgatgacgatgatgatgatgggtttttggaccagccctctcctgatttaaaagcaggatgattctctcgaaggccataataggtcctgcattgccacctggttcattaaa
tctacacagaaaattgacaaggcagagatgatccaggacaacaggtatcatcatccattcaactcaaggtttgaaagttttcattgaaagaataagatgtattgtaacaaactccaaactttatgcatgcagttactacaaccactggactcaactgtaa