

CD300A-CD300a molecule Gene

Reference: PTXBC032352

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD300a molecule; CD300a antigen; CLM-8; CMRF-35-H9; CMRF-35H; CMRF35-H; CMRF35-H9; CMRF35H; CMRF35H9; IGSF12; IRC1; IRC1/IRC2; IRC2; IRp60; CMRF35-like molecule 8; CD300 antigen-like family member A; CMRF35H leukocyte immunoglobulin-like receptor; NK inhi		

For research use only.

GENE DESCRIPTION	CD300a molecule
GENE ID	11314
ORIGIN SPECIES	Human
ACCESSIONS	BC032352
VECTOR	pDONR223
NCBI SYMBOL	CD300A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgctgacctggcctctgttctctctgggtcccagatgttttctctgagcaaatcaggaccctggcggcccccgtgggggatccctgagtgctgagtgccctatgagaaggaaacacagaccctcaacaatactgtgagaccaccacagatttctatgtgacaagatt
gtggagaccaaaagggtcagcaggaaaaaaggaacggcaggtgtccatcagggaacagtcctcgaacctcagcttcacagtgaccctggagaatcacagaggagatgaggcactactgtgtgggtggatcacaccatggctccgagacttcatgatcccgttcgaggttga
gggtctcgtgtcccggcatcaactcaatgacacctgcaagtatcactgcccgaagacctcaacaatcaaaactgcaattccacctgtatcactaccctgttgcagtggtgcccaccacagtgccagatccaggagaaactgaggaaggtggaactcacagctcccctgc
tcctcctcctgctgcatgttctctgttgggtgggctcccctgtagctggaggatgttcagaaatgatcaaaagctgggtgacattcagagctgcccagaacccaagcaggctgccagcagagtgagctgcaactacgaaatcggagctgctgatgtgcctctgcagg
aaaagccagcaccacaaaggaggtgagtggaatacagcactgtggcctccccagggaagaactcactatgcctcgggtgtttgatttaacaccaagatagctgctcagaggcctcgggagaggaaccagattcagattgataaggaagacatag