

TMEM18-transmembrane protein 18 Gene

Reference: PTXBC032379

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 18		

For research use only.

GENE DESCRIPTION	transmembrane protein 18
GENE ID	129787
ORIGIN SPECIES	Human
ACCESSIONS	BC032379
VECTOR	pDONR223
NCBI SYMBOL	TMEM18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctctggatacttcagacggactggactgagccctggctcatgggctggccacctccacgcgctctgctgctcctcactgctgtctctccgaagctacagactacagatcgggcactttctgtgtctagtcatcttagtctactgtgctgaatacatcaatgagcggtgctgcatg
aactggagattatcttgaataaccagatatttcgactccaggggatgttcatttctatagattttcagccccactgctggatgcatgatcattgtgtatgtgggtatggaagacttgaatgtgatgactgacctgaagaatgcacaagagagaagaagaaaaagaaagagaagg
aaagaagactga