



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

LOC440456-similar to pleckstrin homology domain containing, family M (with RUN domain) member 1; adapter protein 162 Gene

Reference: PTXBC032382

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to pleckstrin homology domain containing, family M (with RUN domain) member 1; adapter protein 162
GENE ID	440456
ORIGIN SPECIES	Human
ACCESSIONS	BC032382
VECTOR	pDONR223
NCBI SYMBOL	LOC440456
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctccccgcgaacgggaacgccccactagacaaacgccccagacggcattccagagccgacggcgaggagcagcaggactcgggacagcaaatccccgggggagaagggaagggggaaggcgggacaggggagggccaccctcttggtctccgtacgaccagccagcgcccttcagacgcagcctctctggtccggaccgccggttctctctcgtctcgtcccaatgcccccgaaggctccagacagaggcattcagggtctctctcgtcccccgccagacgcagctcttttgggatgggtgacagcactcttgagcgtcgtctctcttctcgtcgcgtgacgtctcgcgaactctcgcggaccactctcggcaccactcggaaagccccccagggcgatgtctcgggacctaaactctccggccccgccaccgcgccttgggctccctgagggccccagacagactcagggaacgtttgacaaaattgcgaagagagccctacgggagctccgatgttgacacctgctctctcgtccctcagcctcgcgaagcccggtctctcgtcccttggcctgggtcacttagggagggcaggagagatgtgtcagtgtggagcttggagggggccatgcacttcggagaagaccaccgactcatggtttgtgtgagggggccggcgagctccgatt