



Reference: PTXBC032398

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6N1; ATP6N1A; Stv1; VPP1; Vph1; V-type proton ATPase 116 kDa subunit a isoform 1; ATPase, H ⁺ transporting, lysosomal V0 subunit a1; ATPase, H ⁺ transporting, lysosomal non-catalytic accessory protein 1 (110/116kD); H(+)-transporting two-sector ATPase, 1		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal V0 subunit a1
GENE ID	535
ORIGIN SPECIES	Human
ACCESSIONS	BC032398
VECTOR	pDONR223
NCBI SYMBOL	ATP6V0A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr