

MGC40574-hypothetical protein MGC40574 Gene

Reference: PTXBC032412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC40574
GENE ID	285048
ORIGIN SPECIES	Human
ACCESSIONS	BC032412
VECTOR	pDONR223
NCBI SYMBOL	MGC40574
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgctttgccagcggggcgccccgccttcgtccgggtacggaggggacactcccatctagcggcagtgaggctaaagccgctgacagccccactaggaacgatctggagtctgtctagggcactctccctccgcttcacaggttgcacaaaaagcttgggcctgaagctcg
aaggccccggaccgctgggggaagaactag