

NOP16-NOP16 nucleolar protein homolog (yeast) Gene

Reference: PTXBC032424

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOP16 nucleolar protein; NOP16 nucleolar protein homolog; HSPC111; HSPC185; nucleolar protein 16; HBV pre-S2 trans-regulated protein 3; nucleolar protein 16 homolog		

For research use only.

GENE DESCRIPTION	NOP16 nucleolar protein homolog (yeast)
GENE ID	51491
ORIGIN SPECIES	Human
ACCESSIONS	BC032424
VECTOR	pDONR223
NCBI SYMBOL	NOP16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaaggccaaggcgaacccggagcagaagttgttacagtgtcaaccgaagcgtctgaaccggaatgctcgaaggcagcgcgggacgaatgctccacatccgacatgctggaccacgctaatacggtaggcagaacctggccgagatgggttggc
tftggaccccaacaggcggtgcccctccgtaagagaaggtgaagccatggaggtggacatagaggagagcctaagagctgtacggagccctatgtgctgaatgacctggaggcagaagccagcctccagaaaagaagaaatactctgtctgggacctattgactat
gtacgctacatggttagagaaccacggggaggaactataagcagagtggaagacttctcgtattcttgcagaagaggaagatggaggtggagtgactgtttacatcacagctgccccaggtgagggctccccggaccagtgaaagctggagccaggtgttaaggcaaggaggtgc
tgtgtggtccagaggagctggccaggtccatggaatcagaaggttacacacacacgtgcacactccccgctctgggggaaggaaactgttctcagaggtccaatttatattcatctgggggttcacggaaaagccagaacctgtgtttcagggtgggtgatgtaa