



Reference: PTXBC032430

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	integral membrane protein DGCR2/IDD; DGS-C; IDD; LAN; SEZ-12; DiGeorge syndrome critical region protein 2; integral membrane protein deleted in DiGeorge syndrome; DiGeorge syndrome critical region gene 2		

For research use only.

GENE DESCRIPTION	DiGeorge syndrome critical region gene 2
GENE ID	9993
ORIGIN SPECIES	Human
ACCESSIONS	BC032430
VECTOR	pDONR223
NCBI SYMBOL	DGCR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr