

C20orf160-chromosome 20 open reading frame 160 Gene

Reference: PTXBC032455

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20ORF160; C20H20orf160; CCM2-like; dJ310O13.5; cerebral cavernous malformation 2-like; cerebral cavernous malformations 2 protein-like; CCM2 like scaffolding protein		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 160
GENE ID	140706
ORIGIN SPECIES	Human
ACCESSIONS	BC032455
VECTOR	pDONR223
NCBI SYMBOL	C20orf160
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcacgttgcggagtaagctggggccctcgagatccagcagtttgcgatgctgctgctgggagtagccgctgggctgcccacaggactattgcacaggcctgctgaagctctacggagaccggcgcaagttctctcttgggatgcggcccttcacccgaccaggacatc
ggctacttcgagggcttctggaggcgtgggcatccgcgagggcgccatcctcactgacagcttcggccgcatcaagcgcagcatgagctccacgtcggctccgagtgccgagctacgatggcgcgccgagcgcccgaggcacaggccttcaccggctgctgctgacat
cacgcagacatcgagggcgtggccccgatgacgacgacgacgagggatgagccccggggtccaggggcgggagcgacgcccagagacaactacctgtag