

## ZCCHC6-zinc finger, CCHC domain containing 6 Gene

Reference: PTXBC032456

|                |                                                                                                                                                                                                                    |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PAPD6; TUT7; terminal uridylyltransferase 7; PAP associated domain containing 6; TUTase 7; zinc finger CCHC domain-containing protein 6; zinc finger, CCHC domain containing 6; zinc finger CCHC-type containing 6 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger, CCHC domain containing 6                                                            |
| <b>GENE ID</b>          | 79670                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC032456                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | ZCCHC6                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgggagatacagcaaaacatttctggaagcgcactaaagaccgggggactatggatgatgatgacttcagaaggggtcaccaccaagattattataatagatgacctgctaaaggccatggcagtaaaatggaaaaggccctcaaaaaagaagataaccagggaac  
tatgggaataccccagaaaaaggacatgctgtgttccagcaatccatgcatcttaaaacccaatctacagtaaccgcttggaatgacagccacaagatcagagtaagagatggctgtctgatgaacatactggtaactcagacaactggagagaattcaaacctggacctaga  
attcctgttataaacgacaaagaaaagactccttcaagaaaatgaagatgggttaggtgcaagacacaaagagcctgcagaaactgtaagacgactgtttcataaagacctaacaagcctagaaccacgctcagaatggagcaggaagctctgaaaacaagaagcagaggtccaga  
cctagggaagccacggaagactagaatgaggaatgagcaggaatggagactggagggccctgtgacgatgagctgtactttcaacgaaggagctgctaggttacagcagcctgaggagagactgaaagagagactgattgacagggtaaaaaggcgaccacgaaactacct  
acagcaagtagacacctgcagactctgtgatgttttaattgaatccattgcatgtgccataagcatatcaaggaaaagaggcacaagaaaaacattaaggagagaagcaagaggaagagtgctcactacgttaccaccaacacccctccagataaattgcatgtggcattgccattgacaaa  
gtggtacaggaatttgcttacacaatgagaacttgaacagagcgtggaataaacgtatcatggaataatgtgttccaacaaagtaccagattgtccctaaagattatgtggtcatcctgtagcagattgggtttcaaaaatcggatgtaaacattgacatccagtttccagccattatgtct  
cagccagatgctccttactgttcaagaatgtttaagaacagtgactcctttattgatgttgatgcagacttccatgctagggtgccagtggtgtgtgcagagaaaagcaaggctcatttcttaaaactgtttattactaa