

TMC4-transmembrane channel-like 4 Gene

Reference: PTXBC032477

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane channel-like protein 4; testicular tissue protein Li 204; transmembrane channel like 4		

For research use only.

GENE DESCRIPTION	transmembrane channel-like 4
GENE ID	147798
ORIGIN SPECIES	Human
ACCESSIONS	BC032477
VECTOR	pDONR223
NCBI SYMBOL	TMC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagggtggggcggtcagagggatgtggcgctgacaggacgtgtttctgcctcgctccctgggtgctctctctctctctggaatcagatcacttggggggcgactccgaggtgaggactgcaaacctgtggtacaattacaacaactccgtgctgggagactgtc
ctgggcccaggaaatgtacaaactctgctctttgatctgctgactgtcttggcagtcgctgctcatccagtttctagaagaagctcctctggtgctctgctctggggcgctgggtcgtctgctgggggacccaggagttccagggtgccgagagtgctggggtcatctacgcgcagacg
gtggtctgggtggggagttttctgccccttactgcccctgttaacacgggtcaagtctcgtctgtttctacctgaagaagcttaacctcttctccacctgtccccgggtgcccgcaccttccgggctccgcccgaatttcttttcccttgggtctctctgggtctggtccatctccagct
tccccgtttacagcatcttctgatcccgcttctaaagcttgggtccattccggggcagtcgctcatctgggccagatccctgagctatcttccagctccctgagaccaccagaatttctcttcttctggggaccaggttttctgctgcccctctgctgatctccagatctgatg
gcgtacactgtggtctggttaactctacggacgctcatctgagctcaaacgtcagagacagcggagggcgcaataaagtcttctggtgacggcgctgtggtgctgacctccaccaaacggctcttga