



Reference: PTXBC032478

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAPKK3; MEK3; PRKMK3; SAPKK-2; SAPKK2; dual specificity mitogen-activated protein kinase kinase 3; MAP kinase kinase 3; MAPK/ERK kinase 3; MAPKK 3; MEK 3; SAPK kinase 2; stress-activated protein kinase kinase 2; mitogen-activated protein kinase kinase 3		

For research use only.

GENE DESCRIPTION	mitogen-activated protein kinase kinase 3
GENE ID	5606
ORIGIN SPECIES	Human
ACCESSIONS	BC032478
VECTOR	pDONR223
NCBI SYMBOL	MAP2K3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1