

SPIN3-spindlin family, member 3 Gene

Reference: PTXBC032490

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPIN-3; TDRD27; bA445O16.1; spindlin-3; bA445O16.1 (DXF34); spindlin-like protein 3; spindlin family member 3		

For research use only.

GENE DESCRIPTION	spindlin family, member 3
GENE ID	169981
ORIGIN SPECIES	Human
ACCESSIONS	BC032490
VECTOR	pDONR223
NCBI SYMBOL	SPIN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagacccgtttgaaaggcagctgcaggcagcggtccaggacggcgctgccacggcagtggtgtgttaccatgataaagggaaggctgcacacaagaagcataggagccgacccacctccagcctcgggggaacatcgtgggctgcagaattcagcagcgatgga
aagatggagatgaaccttaacacagtggaagggaaccgttctggatcagctccttga