



Reference: PTXBC032503

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APC2; anaphase-promoting complex subunit 2; cyclosome subunit 2; anaphase promoting complex subunit 2		

For research use only.

GENE DESCRIPTION	anaphase promoting complex subunit 2
GENE ID	29882
ORIGIN SPECIES	Human
ACCESSIONS	BC032503
VECTOR	pDONR223
NCBI SYMBOL	ANAPC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]