



Reference: PTXBC032504

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HIPER1; MINPP2; MIPP; multiple inositol polyphosphate phosphatase 1; 2,3-BPG phosphatase; 2,3-bisphosphoglycerate 3-phosphatase; inositol (1,3,4,5)-tetrakisphosphate 3-phosphatase; ins(1,3,4,5)P(4) 3-phosphatase; multiple inositol polyphosphate histidine		

GENE DESCRIPTION	multiple inositol polyphosphate histidine phosphatase, 1
GENE ID	9562
ORIGIN SPECIES	Human
ACCESSIONS	BC032504
VECTOR	pDONR223
NCBI SYMBOL	MINPP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]