

NGB-neuroglobin Gene

Reference: PTXBC032509

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	neuroglobin
GENE ID	58157
ORIGIN SPECIES	Human
ACCESSIONS	BC032509
VECTOR	pDONR223
NCBI SYMBOL	NGB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcggccggagccgagctgacccggcagagctggcgggagtgagccgagccgctggagcacggcaccgtctgtttgccaggctgttggccctggagcctgacgtgctgccctcttcagtaacaactgccgagcttctcagccagaggagctgtctctcctgcctga
gttcctggaccacatcaggaaggtgatgctctgtgattgatgctgacgaatgtggaagacctgtcctcactggaggagtagcttgcagccctggcagggaagcaccggcagtggtgtggaagctcagctccttcgacagtggtgagctctctacatgctggagaaggtgc
tgggcccctgccttcacaccagccacacgggctgcctggagccaactctacggggccgtagtgaggccatgagtcgagcctgggatggcgagtaa