

FAM86C-family with sequence similarity 86, member C Gene

Reference: PTXBC032519

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM86C; protein FAM86C1; family with sequence similarity 86, member C; family with sequence similarity 86 member C1		

For research use only.

GENE DESCRIPTION	family with sequence similarity 86, member C
GENE ID	55199
ORIGIN SPECIES	Human
ACCESSIONS	BC032519
VECTOR	pDONR223
NCBI SYMBOL	FAM86C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggaggagaacgcggggagcgaactctgtcagagtttcaagcggcgttctctgagcgcgcggcctgcgtcttccgctggcagagcttagaagcaaaagtaagagactcatcagattctgagctgctgcgggatatttcagaagacgtgctgtattgccagaagccatcgtgtcgtggtcgggtcctcgggagcgtgctgctgcccggagcaccagcgggctcctcaattctacatggcccttaccgtctgcaaccagagatgtgccagctgttcaccaccgagctatgctggactgggagcagatgggaagcgggaagctcatgaccagaaactgttcctacagagagcacttgagatggcaagctga