

DTWD1-DTW domain containing 1 Gene

Reference: PTXBC032535

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDS009; DTW domain-containing protein 1; x 009 protein; DTW domain containing 1		
For research use only.			

GENE DESCRIPTION	DTW domain containing 1
GENE ID	56986
ORIGIN SPECIES	Human
ACCESSIONS	BC032535
VECTOR	pDONR223
NCBI SYMBOL	DTWD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcctccacgtatatttcaaacgaagtgaagaaatgttcaaaattgtggaacaaacagtcacaaactctccatagcttcagaagatcccttcaaaacttatgtttagcatctcaagaagttcttcaaaagctcagcaaaagtgaggatcaaaatgtctcaaatgtgtgtgtgttccagaatgttctactgctatactgttatgttccagttgaaaatgtacattgaacagattccattgtgaagcttcattgaagattgacatcattaaacatccaaatgaacagatggcaaaagttctctatactgcaaaactcttagcactgaatttgaacatttacagtatccgtgtattccagaatgaagaaaggacatgaagttgcaactcattttctggacctcagctctatctcaataaaagatattttttcatctgcaaaaaaggattcaaaataatgttagaggcaaaatgatgacctgacaagccatcttttaaacgcaaaagaaactgaagaacaagagttctgtgattgaaatgacagcaagtgcaaggcacaacactgaaaaaattattttatagatgacctggaacaaacaacaaatattcactgatgagcgacttcaagggtgtttacaaagttgagttgaaaacaagaaaaactgttttggcgccatcaaaaggaaagccagatactttcttctacaattgaagccatttactactttctgttagactaccatactgatatattaaagagaatacagaggccaatgacaatcttttttttctattctttgtaccagttgataaagaatgccaaatgctctggagataaggaaacaggaaacttacacattag