

LOC729603-calcium binding protein P22 pseudogene Gene

Reference: PTXBC032536

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	calcium binding protein P22 pseudogene
GENE ID	729603
ORIGIN SPECIES	Human
ACCESSIONS	BC032536
VECTOR	pDONR223
NCBI SYMBOL	LOC729603
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtctcaggcctccacttactgcacgatgaagagttttaggagatcaagaaggagactggctttccacagtcacacgtctgtacagccggttcagcaacctggacaaaggagagacaggacgattccaggggattccagaacttccatcaaccagtgaggact
ggatcatcaatgccttccctccagaggagagaccaggttaa