

POU2AF1-POU class 2 associating factor 1 Gene

Reference: PTXBC032549

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BOB1; OBF-1; OBF1; OCAB; POU domain class 2-associating factor 1; B-cell-specific coactivator OBF-1; BOB-1; OCA-B; OCT-binding factor 1; POU class 2 associating factor 1		

For research use only.

GENE DESCRIPTION	POU class 2 associating factor 1
GENE ID	5450
ORIGIN SPECIES	Human
ACCESSIONS	BC032549
VECTOR	pDONR223
NCBI SYMBOL	POU2AF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctggcaaaaacccacagctccggagcaagccccagccccggccggccataccaggcgctcgtgtgaaaggagccagtgaaggaaactgctgaggagggaagcaggccacgcagcagctggggcagcacctgcacctagggcggtgctgccccatcagccccctggc
gacctacaccacagtggtctctctcctggacatggaaggttctgtgtcagtgacagaggagggctccctgtgtgcccgtggtctctccagccccacccggccacctgcagccccctggccccatggacaccttacaccgagtagtgcacctgaagctgtcagctgccccctc
ctcagctgacatgtatgtgcagccccgtgtgccccagctacacgggtgtggggccctcctcagtggtgacctatgcctctcggccactcatcaacaatgtcacgacaagaagctccgccacggccagctggggccccctgctggaggggccagagcaccaggcaccctcactatctcc
gtggcctcagcccccttccacactaccacctccacctgcagtagcagcctccggccccagccctacctggccccagttgtccagctcccatctctatccagagccagctccttcagacatggaagacccagagaagccgagcagctcgttgaccatcgacaagctgttttgag
gaagaggatagcagcctatgcgttaaccacactctctgtggaaggttttag