

ZNF702P-zinc finger protein 702 pseudogene Gene

Reference: PTXBC032590

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	zinc finger protein 702 pseudogene
GENE ID	79986
ORIGIN SPECIES	Human
ACCESSIONS	BC032590
VECTOR	pDONR223
NCBI SYMBOL	ZNF702P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagaaaaatcttccaatgtaatgagagtggaagcctttaattgtagctcactgttaaaaaaatgtcagataatccatttaggagagaaaaataataatgtgatataatgtgcaaggtctttaatcagaagcgataccttgcataccatcatagatgtcacactggtgagaacaccttacaagtgtaatcagtggtgcaagaccttcagttacaagtcaccccttgaattcacaaggcaattcactggagagaaacctcacaagtgaatgaatgtgcaaggttttaatacaaaagcatatcttgcagtcacatagacttcatactggagagaaccttacaagtgtgaagaatgtgacaaagtgttagctgctaa