

IFI27L2-interferon, alpha-inducible protein 27-like 2 Gene

Reference: PTXBC032626

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM14A; ISG12B; TLH29; interferon alpha-inducible protein 27-like protein 2; ISG12(b) protein; family with sequence similarity 14, member A; interferon-stimulated gene 12b protein; pIFI27-like protein; interferon alpha inducible protein 27 like 2		

For research use only.

GENE DESCRIPTION	interferon, alpha-inducible protein 27-like 2
GENE ID	83982
ORIGIN SPECIES	Human
ACCESSIONS	BC032626
VECTOR	pDONR223
NCBI SYMBOL	IFI27L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgaaacgggcagctgctgctgcagtgaggagccctggcagtgaggcctgtgcccgtggtgctcagtgccatgggcttactggggcaggaatcgcgcgtctccatagcagccaagatgatgtccgcagcagccattgccaacgggggtggtgttctgctggggagcctg
gtggctactctgcagtcctgggggcagctggactctccacatcatccaatcctctgacctctgttgggtcagtggtggggcctgcttgggggaattcacccttcttctctccagctgaacccgaggctaaagaagatgaggcaagagaaaatgtacccaaggtgaacctcaaaa
ccccactcaagtcagagaacatgaggaataa