



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF41-ring finger protein 41 Gene

Reference: PTXBC032637

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

FLRF; NRDP1; SBBI03; E3 ubiquitin-protein ligase NRDP1; fetal liver ring finger; neuregulin receptor degradation protein-1; ring finger protein 41, E3 ubiquitin protein ligase; ring finger protein 41

For research use only.

GENE DESCRIPTION	ring finger protein 41
GENE ID	10193
ORIGIN SPECIES	Human
ACCESSIONS	BC032637
VECTOR	pDONR223
NCBI SYMBOL	RNF41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggatgatgtaaacccgtttccagggggatggtgacgaagaatctatctgccctatttgcagttggagctctggaggagccagtagcacggcactcattggaacatctttctcaacgctgcatcaccacagtggttctctcagcaacagacatgtccagtgggaccgtagtggtgacggt
gcgccatctgcgccagtgccggatcatcggaacatgtgtgaacaagtctgacgaacgtgtgttcggctgtagtccgtttgcggcttgaccaacatcatgtcaacatcagcagatgtgacgacaaccccgccgctgtgacctgtgaacacggcgctgtgacacggtgtgacacggt
agatgcaccaagaatgagctgcgcccaataaactgcattaaagcactgcagttgtgtacagcagcagacagcagtcgacagatctggagaagacctgcagctgcaaaccaacacagctgcggcagcaagaagcagatccagctgtcaatgaacatgctgtgcaatccgc
agtgtcaaccccaaccttcagaacctggaggagacaattgaatacaacgagatcctagatgggtggaactccctccagccagcaagatgtagcccgctggggagagctgacatcgtactcctgatctgtctgccacggctgtaatacaagcgtccctgtgtggagagtggtctgtcctcttat
tgtcaacgagctgattgaaatgccccagcagctagctgtgccccaggctgtgtgccacatagacatgacagatgacacgagcgtactatgagaacatctgtgcccaagcgcataccctgtcgaacagcgtctgttctgtgatggcctgtgagaaccagacatggggatgacatgggt
caagaagccagcctcttctcatgatatatttcccatgcctctggaagagataata

More info: Please contact us for any question: contact@proteogenix.fr