



Reference: PTXBC032674

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBH1; Fbx18; hFBH1; F-box DNA helicase 1; F-box only protein, helicase, 18; F-box protein, helicase, 18		

For research use only.

GENE DESCRIPTION	F-box protein, helicase, 18
GENE ID	84893
ORIGIN SPECIES	Human
ACCESSIONS	BC032674
VECTOR	pDONR223
NCBI SYMBOL	FBXO18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]