



Reference: PTXBC032688

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGTRL1; APJR; HG11; apelin receptor; APJ (apelin) receptor; APJ receptor; G protein-coupled receptor APJ; G-protein coupled receptor HG11; HG11 orphan receptor; angiotensin II receptor-like 1; angiotensin receptor-like 1		

For research use only.

GENE DESCRIPTION	apelin receptor
GENE ID	187
ORIGIN SPECIES	Human
ACCESSIONS	BC032688
VECTOR	pDONR223
NCBI SYMBOL	APLNR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]