

GPR55-G protein-coupled receptor 55 Gene

Reference: PTXBC032694

Product type: DNA Storage Condition: -20°C
Size/Quantity: 2 ug Delivery Condition: +4°C, Blue Ice
Synonyms: LPIR1; G-protein coupled receptor 55; G protein-coupled receptor 55
For research use only.

GENE DESCRIPTION	G protein-coupled receptor 55
GENE ID	9290
ORIGIN SPECIES	Human
ACCESSIONS	BC032694
VECTOR	pDONR223
NCBI SYMBOL	GPR55
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcagcaaaaccagtgaggactgctgtttgacggtgtcaacgagctgatgaaaaccctacagttgacgtccacatcccacctctgctggcctgctcctcaacctgctggccatccatgcttcagcaccctcctaagaacagtgcccattatgctccacctccatc
acatgatcaacctggcagctctttgacctgctgctggctcctcccattcaagatggtctgtcccaggtacagtcocccctcccgctcctgtgacacctgggtgagtgctttacttctgctcagcatgtacggaagcgtcttccacatctgcttcacagcatggaccggtcttgccatccgt
taccgctactggtgagccactccggtccccaggaagatctttggatctgctgaccatctgggtcctggtgaccggaagcaccctctacagtttccatgggaaagtggaaaaatcatgtctccacaacatgctgatatacctggagcgcgaaggtcttctcccgctgg
agggtttggcttctcctccatgggcatcatggcttctgctgctccaggagcaccacatctgctgggcccggagaccacaccaggaactgggtgcagcagaagcctgcatctacagcatgcagccagcctggctgctgctgctcctcccaagtcaccctgggttct
cctgcagttcctggtgagaacagctttagtagtgcagagccaagcagagcatcagcttcttctgcaattgtccatgtgttttccaacgtcaactgctgctggatgtttctgctactactttgcatcaagaattcgcgatgaacatcagggccaccggcctccaggtccagctg
gtctcgaggacaccagatctcccgggctaa