

RPL30-ribosomal protein L30 Gene

Reference: PTXBC032700

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L30; 60S ribosomal protein L30; ribosomal protein L30		

For research use only.

GENE DESCRIPTION	ribosomal protein L30
GENE ID	6156
ORIGIN SPECIES	Human
ACCESSIONS	BC032700
VECTOR	pDONR223
NCBI SYMBOL	RPL30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggccgcaaagaagacgaaaaagtcgtggagtcgatactctaggtccaactcgttatgaaaagtgggaagtagctcctgggtacaagcagactctgaagatgatcagacaaggcaaaagcgaattggtcattctcgctaacaactgccagctttgaggaaatctgaaatag
agtactatgctatgttggctaaaactggtgtccatcactacagtggcaataatattgaactgggcacagcatgcggaaaatactacagagtgtgcacactggctatcattgatccaggtgactctgacatcattagaagcatgccagaacagactggtgaaaagtaa