



Reference: PTXBC032743

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone acetyltransferase KAT2A; GCN5; GCN5L2; PCAF-b; hGCN5; GCN5 (general control of amino-acid synthesis, yeast, homolog)-like 2; General control of amino acid synthesis, yeast, homolog-like 2; K(lysine) acetyltransferase 2A; STAF97; general control of		

For research use only.

GENE DESCRIPTION	K(lysine) acetyltransferase 2A
GENE ID	2648
ORIGIN SPECIES	Human
ACCESSIONS	BC032743
VECTOR	pDONR223
NCBI SYMBOL	KAT2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]