



Reference: PTXBC032778

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAO; peroxisomal N(1)-acetyl-spermine/spermidine oxidase; peroxisomal N1-acetyl-spermine/spermidine oxidase; polyamine oxidase (exo-N4-amino); polyamine oxidase		

For research use only.

GENE DESCRIPTION	polyamine oxidase (exo-N4-amino)
GENE ID	196743
ORIGIN SPECIES	Human
ACCESSIONS	BC032778
VECTOR	pDONR223
NCBI SYMBOL	PAOX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr