

TFF2-trefoil factor 2 Gene

Reference: PTXBC032820

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SML1; trefoil factor 2; spasmodic; spasmodic polypeptide; spasmodic protein 1		
For research use only.			

GENE DESCRIPTION	trefoil factor 2
GENE ID	7032
ORIGIN SPECIES	Human
ACCESSIONS	BC032820
VECTOR	pDONR223
NCBI SYMBOL	TFF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacggcgagagcggcagctcctggcagcgtcctcctggggtatgtgcccggcggggagtgagaaaccctccccctgagtgctcaggtgagccccataacaggacgaactcgggtccctggaatcaccagtgaccagtggttgacaatggatgctgttcgac
tccagtgctcactgggtccctgggtgtttccaccccccccaagcagagtcggatcagtgctcagtgaggtcagaccgaagaaactgtggctacccgggcacagccccgagggaatgcgcctctcggaagtgccttccaactcatttgaagtgcctgtcttccg
aagctgtggaagactgccattactaa