



Reference: PTXBC032856

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TTC18; cilia- and flagella-associated protein 70; TPR repeat protein 18; tetratricopeptide repeat domain 18; tetratricopeptide repeat protein 18; tetratricopeptide repeat-containing protein (LOC118491); cilia and flagella associated protein 70		

For research use only.

GENE DESCRIPTION	tetratricopeptide repeat domain 18
GENE ID	118491
ORIGIN SPECIES	Human
ACCESSIONS	BC032856
VECTOR	pDONR223
NCBI SYMBOL	TTC18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]