



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF345-zinc finger protein 345 Gene

Reference: PTXBC032863

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HZF10; zinc finger protein 345; zinc finger protein HZF10

For research use only.

GENE DESCRIPTION	zinc finger protein 345
GENE ID	25850
ORIGIN SPECIES	Human
ACCESSIONS	BC032863
VECTOR	pDONR223
NCBI SYMBOL	ZNF345
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr