

AQP5-aquaporin 5 Gene

Reference: PTXBC032946

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AQP-5; PPKB; aquaporin-5; aquaporin 5		
For research use only.			

GENE DESCRIPTION	aquaporin 5
GENE ID	362
ORIGIN SPECIES	Human
ACCESSIONS	BC032946
VECTOR	pDONR223
NCBI SYMBOL	AQP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaaggaggtgtgtctccgtggccttctcaaggccgtgttcgagagtttggccaccctcatctctcttcttggcctgggctcggccctcaagtgccgtcggcgctgcctaccatcctgcagatcgcgtggcggttggcctggccataggcacgtggccaggccctggg
acccgtgagcggcggccacatcaaccccgccatcccttgggtgggcaaccagatctcgtctcggccttctctctacgtggcggccagctgggtggcgccattgccgggctggcatcctctacgtgtggcaccgtcaatgccggggaatctggcgtcaacg
cgctcaacaacaacaacgagggccagggccatgggtgggtgagctgattctgacctccagctggcactctgcatcttcgctccactgaetcccgccgacacagccctgtgggtccccagccctgtccattggcctgtctgtcaccctggccaccttgcggaatctacttactggct
gtcccatgaaccagcccgctcttttggcctgcgggtggtcatgaatcggttcagccccgctcactgggttttctgggtaggggccatcgtggggcggtcctggtgccatcttacttctactctcttccccaaactccctgagcctgagtgagcgtgtggccatcatcaaggcacgtat
gagcctgacgaggactgggagggagcagcgggaagagcgggaagaagaccatggagctgaccaccgctga