



Reference: PTXBC032953

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	fibrinogen C domain-containing protein 1; fibrinogen C domain containing 1		
For research use only.			

GENE DESCRIPTION	fibrinogen C domain containing 1
GENE ID	84929
ORIGIN SPECIES	Human
ACCESSIONS	BC032953
VECTOR	pDONR223
NCBI SYMBOL	FIBCD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]