

C1orf111-chromosome 1 open reading frame 111 Gene

Reference: PTXBC032957

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf111; C1orf111; HSD20; spermatogenesis associated 46		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 111
GENE ID	284680
ORIGIN SPECIES	Human
ACCESSIONS	BC032957
VECTOR	pDONR223
NCBI SYMBOL	C1orf111
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaactctcactctcagcatctctggacctccaatctctcctcctcctgagtgctttcccgacattatgtctctgtgccaccagcctgccagacattgcaaagacagcagtagccactgaggcatccagccagctcagggcctgccacccagtagcaaaagcatcattgtcag
gcaaggatagacagaacacagtgctctcaccagactgcagcttgggggacacccagcacggagagaagctgaggcggaaactgactatctaccggcctgtctctccctacagctactctgtgtgtgagacaaagagagccagctggaggcctatgacttccagaggtgcagcag
gatgaggcgaagtgggacaaactgcttctgaggacatggctgagaacatctgtctcctctctcctcccgagaaacacttgcctctgagaagccaccaagaatccaggcatggcctggactccatcacatcccaggacatcctaattggctccagatggcaccagcacagcagaatg
gctacaagtgcgtggcctgtgcccgaatgtacccaccctggacttctcctcaagagccacatcaagaggggcttcaggagggttcagctgcaaggtgtactaccacaagctcaaaagcctctggagcaaggagcagaaggcccgctgggagacaggtctctccggcagctgcc
aggcctcaatagctcgtgtaacaccttaggcaaatggcgtgaagcctacttatgtctctag