

RP4-662A9.2-hypothetical protein MGC34034 Gene

Reference: PTXBC032958

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein MGC34034		
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC34034
GENE ID	154089
ORIGIN SPECIES	Human
ACCESSIONS	BC032958
VECTOR	pDONR223
NCBI SYMBOL	RP4-662A9.2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagtggcttcaagagctcttcagcaaccaggcaccttcagggtcttcttcccaattctagtgctagtaactcttcctcatcatccaacatgggtaccagagctccagacatcacgtgtacatctgagttccaggaaaccagacagagggagaaaaagggcatgcacgtcactt
aaagaagacttctaagaagaaccacctaataatcccccaactagtcctcaagctcacttggctgctttcaagcctggaataatgcagctcttaagctgaaaaatgctatgacaaatcagagattgatcatgctacaatgacaactggaaaacagaaatgaaccatataagccaagcactca
cgcaatcttctcctccagttgacttaagaacacttgaagaaaatgtggcgctttag