



Reference: PTXBC033068

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-BETA13; protocadherin beta-13; protocadherin beta 13		
For research use only.			

GENE DESCRIPTION	protocadherin beta 13
GENE ID	56123
ORIGIN SPECIES	Human
ACCESSIONS	BC033068
VECTOR	pDONR223
NCBI SYMBOL	PCDHB13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr